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---

## AIRONE 200

---

**BIOQUIMICA**

BioTécnica Ind e Com. Ltda.

Tel / Fax: +55 35 3214-4646 Varginha MG Brasil.

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Temperatura: 37 °C

TEST NAME	ALB
DECIMAL	1
MEASURE UNIT	g/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	3
REAGENT VOLUME	
REAGENT 1	600
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	6.0
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	AMI
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	1038
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	6950

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TEST NAME	ALT
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	40
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	350
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	-1746

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TEST NAME	AST
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	40
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	350
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	-1746

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TEST NAME	BDI
DECIMAL	2
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	6
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	15
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFERENTIAL	YES
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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TEST NAME	BTO
DECIMAL	2
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	6
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	15
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFERENTIAL	YES
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	CA
DECIMAL	2
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	20
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	CA
DECIMAL	2
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	20
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	578
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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TEST NAME	CK
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	16
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	1000
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	8095

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TEST NAME	CK MB
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	16
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	600
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	15
MEASURE	120
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	6750

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Temperatura: 37 °C

TEST NAME	CLO
DECIMAL	0
MEASURE UNIT	mEq/L
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	3
REAGENT VOLUME	
REAGENT 1	600
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	120
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	500
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	COL
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	800
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	CREA
DECIMAL	2
MEASURE UNIT	mg/dL
REACTION TYPE	FIXED TIME
SAMPLE VOLUME (µL)	40
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	12
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	30
INCUBATION 2	0
LAG PHASE	15
MEASURE	60
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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TEST NAME	LDH
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	4000
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	-8095

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Temperatura: 37 °C

TEST NAME	FER
DECIMAL	0
MEASURE UNIT	µg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	16
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	500
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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TEST NAME	FAL
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
<b>REAGENT VOLUME</b>	
REAGENT 1	400
REAGENT 2	0
<b>ABSORBANCE RANGE (mABS)</b>	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	700
CONTAMINATING	YES
REAGENT BLANKING	NO
DIFERENTIAL	NO
<b>FILTERS (nm)</b>	
FILTER 1	405
FILTER 2	NONE
<b>TIMES (SEC)</b>	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
<b>MEASUREMENT TYPE</b>	
CALIBRATE	FACTOR
FACTOR	2757

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Temperatura: 37 °C

TEST NAME	FUV
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	15
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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TEST NAME	GGT
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	40
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	1000
CONTAMINATING	YES
REAGENT BLANKING	NO
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	1158

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Temperatura: 37 °C

<b>TEST NAME</b>	GLI
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
<b>REAGENT VOLUME</b>	
REAGENT 1	400
REAGENT 2	0
<b>ABSORBANCE RANGE (mABS)</b>	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	400
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
<b>FILTERS (nm)</b>	
FILTER 1	510
FILTER 2	NONE
<b>TIMES (SEC)</b>	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	2
<b>MEASUREMENT TYPE</b>	
CALIBRATE	STANDARD
FACTOR	0

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TEST NAME	HDL
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	40
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	150
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	HDL
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	FIXED TIME
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	300
REAGENT 2	100
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	180
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	600
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.20
MIX 2	0.20
INCUBATION 1	300
INCUBATION 2	5
LAG PHASE	3
MEASURE	120
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	LDL
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	FIXED TIME
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	300
REAGENT 2	100
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	400
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	600
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.20
MIX 2	0.20
INCUBATION 1	300
INCUBATION 2	5
LAG PHASE	3
MEASURE	120
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	MAG
DECIMAL	2
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	4.5
CONTAMINATING	NO
REAGENT BLANKING	
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	PRO
DECIMAL	2
MEASURE UNIT	g/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	12
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	PRO-U
DECIMAL	2
MEASURE UNIT	mg/L
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	8
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	700
LINEARITY LIMIT	4000
CONTAMINATING	NO
REAGENT BLANKING	
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	TGR
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	1000
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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TEST NAME	URE
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	3
<b>REAGENT VOLUME</b>	
REAGENT 1	300
REAGENT 2	300
<b>ABSORBANCE RANGE (mABS)</b>	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	200
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFERENTIAL	NO
<b>FILTERS (nm)</b>	
FILTER 1	620
FILTER 2	NONE
<b>TIMES (SEC)</b>	
MIX 1	0.00
MIX 2	0.20
INCUBATION 1	300
INCUBATION 2	300
LAG PHASE	3
MEASURE	1
<b>MEASUREMENT TYPE</b>	
CALIBRATE	STANDARD
FACTOR	0

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TEST NAME	UUV
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0.800
MAXIMUM	2000
LINEARITY LIMIT	250
CONTAMINANTING	NO
REAGENT BLANKING	NO
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	30
INCUBATION 2	0
LAG PHASE	10
MEASURE	60
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

Temperatura: 37 °C