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AIRONE 200

TURBIDIMETRIA

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Temperatura: 37 °C

TEST NAME	AGP
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (μ L)	20
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	200
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	APO-A
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	8
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	300
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	APO-B
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	10
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	300
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	ASO
DECIMAL	1
MEASURE UNIT	UI/mL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	200
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.2
MIX 2	0.00
INCUBATION 1	180
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	

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Temperatura: 37 °C

TEST NAME	C3
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	8
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	500
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	C4
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	20
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	100
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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TEST NAME	FR
DECIMAL	2
MEASURE UNIT	UI/mL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	100
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.2
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	

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TEST NAME	FR
DECIMAL	2
MEASURE UNIT	UI/mL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	100
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.2
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	MULTI STD
FACTOR	

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Temperatura: 37 °C

TEST NAME	FER
DECIMAL	1
MEASURE UNIT	µg/L
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	40
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	300
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.2
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	HBA1C
DECIMAL	1
MEASURE UNIT	%
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	8
REAGENT VOLUME	
REAGENT 1	300
REAGENT 2	100
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	1000
LINEARITY LIMIT	16
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.2
MIX 2	0.2
INCUBATION 1	300
INCUBATION 2	300
LAG PHASE	3
MEASURE	2
MEASUREMENT TYPE	
CALIBRATE	CUBIC
FACTOR	0

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Temperatura: 37 °C

TEST NAME	IGA
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	600
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	IGG
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	3000
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	IGM
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	20
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	500
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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Temperatura: 37 °C

TEST NAME	MAL
DECIMAL	1
MEASURE UNIT	mg/L
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	3
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	80
CONTAMINANTING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.2
MIX 2	0.00
INCUBATION 1	180
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	

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Temperatura: 37 °C

TEST NAME	PCR
DECIMAL	2
MEASURE UNIT	mg/L
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	4
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	150
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.2
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	

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Temperatura: 37 °C

TEST NAME	TFR
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	End Point
SAMPLE VOLUME (µL)	8
REAGENT VOLUME	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	500
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0