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COBAS MIRA S / PLUS

TURBIDIMETRIA

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GENERAL		CONVERSION FACTOR	
MEASUREMENT MODE	ABSORB	OFFSET	1.000
REACTION MODE	D-R-S	TEST RANGE LOW	0.000
CALIBRATION MODE	SLOP AVG	HIGH	30 mg/dL
REAGENT BLANK	REAG/DIL	NORM RANGE LOW	200 mg/dL
CLEANER	NO	HIGH	40 mg/dL
WAVELENGTH	340 nm	NUMBER OF STEPS	150 mg/dL
DECIMAL POSITION	1	CALCULATION STEPS A	1
UNIT	mg/dL	READING FIRST	1
		LAST	25
ANALYSIS		CALIBRATION	
DILUTION NAME	NaCl	CALIB. INTERVAL	ON REQUEST
FACTOR	10	TIME	NO
TIME	NO	BLANK RANGE LOW	NO
STD	MAIN INDIRECT	HIGH	NO
MAIN STD	##	REAGENT RANGE LOW	NO
POST DIL. FACTOR	2	HIGH	NO
CONC. FACTOR	NO	STANDARD POS	*
SAMPLE CYCLE	1	STD-1	@
VOLUME	12,5 µL	STD-2	
DILUTION NAME	NaCl	STD-3	
VOLUME	00 µL	REPLICATE	DUPL
REAGENT CYCLE	1	DEVIATION	15%
VOLUME	250 µL	CONTROL	
CALCULATION		SAMPLE LIMIT	NO
SAMPLE LIMIT	NO	CS1 POS	*
REACTION DIRECTION	INCREASE	CS2 POS	*
CHECK	ON	CS3 POS	*

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GENERAL		CONVERSION FACTOR	
MEASUREMENT MODE	ABSORB	OFFSET	1.000
REACTION MODE	D-R-S	TEST RANGE LOW	0.000
CALIBRATION MODE	SLOP AVG	HIGH	10 mg/dL
REAGENT BLANK	REAG/DIL	NORM RANGE LOW	300 mg/dL
CLEANER	NO	HIGH	110 mg/dL
WAVELENGTH	340 nm	NUMBER OF STEPS	210 mg/dL
DECIMAL POSITION	1	CALCULATION STEPS A	1
UNIT	mg/dL	READING FIRST	ENDPOINT
		LAST	1
			25
ANALYSIS		CALIBRATION	
DILUTION NAME	NaCl	CALIB. INTERVAL	ON REQUEST
FACTOR	10	TIME	NO
TIME	NO	BLANK RANGE LOW	NO
STD	MAIN INDIRECT	HIGH	NO
MAIN STD	##	REAGENT RANGE LOW	NO
POST DIL. FACTOR	2	HIGH	NO
CONC. FACTOR	NO	STANDARD POS	*
SAMPLE CYCLE	1	STD-1	@
VOLUME	5,0 µL	STD-2	
DILUTION NAME	NaCl	STD-3	
VOLUME	00 µL	REPLICATE	DUPL
REAGENT CYCLE	1	DEVIATION	15%
VOLUME	250 µL		
CALCULATION		CONTROL	
SAMPLE LIMIT	NO	CS1 POS	*
REACTION DIRECTION	INCREASE	CS2 POS	*
CHECK	ON	CS3 POS	*

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GENERAL		CONVERSION FACTOR	
MEASUREMENT MODE	ABSORB	OFFSET	1.000
REACTION MODE	D-R-S	TEST RANGE LOW	0.000
CALIBRATION MODE	SLOP AVG	HIGH	10 mg/dL
REAGENT BLANK	REAG/DIL	NORM RANGE LOW	300 mg/dL
CLEANER	NO	HIGH	60 mg/dL
WAVELENGTH	340 nm	NUMBER OF STEPS	155 mg/dL
DECIMAL POSITION	1	CALCULATION STEPS A	1
UNIT	mg/dL	READING FIRST	ENDPOINT
ANALYSIS		LAST	25
DILUTION NAME	NaCl	CALIBRATION	
FACTOR	10	CALIB. INTERVAL	ON REQUEST
TIME	NO	TIME	NO
STD	MAIN	BLANK RANGE LOW	NO
MAIN STD	##	HIGH	NO
POST DIL. FACTOR	2	REAGENT RANGE LOW	NO
CONC. FACTOR	NO	HIGH	NO
SAMPLE CYCLE	1	STANDARD POS	*
VOLUME	6,0 µL	STD-1	@
DILUTION NAME	NaCl	STD-2	
VOLUME	00 µL	STD-3	
REAGENT CYCLE	1	REPLICATE	DUPL
VOLUME	240 µL	DEVIATION	15%
CALCULATION		CONTROL	
SAMPLE LIMIT	NO	CS1 POS	*
REACTION DIRECTION	INCREASE	CS2 POS	*
CHECK	ON	CS3 POS	*

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GENERAL		TEST RANGE	LOW	0 IU/mL
MEASUREMENT MODE	ABSORB		HIGH	800 IU/mL
REACTION MODE	R-S	NORM RANGE	LOW	0 IU/mL
CALIBRATION MODE	SLOP AVG		HIGH	200 IU/mL
REAGENT BLANK	REAG/DIL	NUMBER OF STEPS		1
CLEANER	NO	CALCULATION STEPS A		ENDPOINT
WAVELENGTH	550 nm	READING FIRST		1
DECIMAL POSITION	1	LAST		6
UNIT	IU/mL			
ANALYSIS		CALIBRATION		
POST DIL. FACTOR	2	CALIB. INTERVAL		ON REQUEST
CONC. FACTOR	NO	TIME		NO
SAMPLE CYCLE	1	BLANK RANGE	LOW	NO
VOLUME	3,0 µL		HIGH	NO
DILUTION NAME	H2O	REAGENT RANGE	LOW	NO
VOLUME	20.0 µL		HIGH	NO
REAGENT CYCLE	1	STANDARD POS		*
VOLUME	300 µL	STD-1		@
		STD-2		
		STD-3		
CALCULATION				
SAMPLE LIMIT	NO			
REACTION DIRECTION	INCREASE	REPLICATE		DUPL
CHECK	ON	DEVIATION		15%
CONVERSION FACTOR	1.000	CS1 POS		*
OFFSET	0.000	CS2 POS		*
		CS3 POS		*

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GENERAL		TEST RANGE	LOW	0 IU/mL
MEASUREMENT MODE	ABSORB		HIGH	800 IU/mL
REACTION MODE	R-S	NORM RANGE	LOW	0 IU/mL
CALIBRATION MODE	SLOP AVG		HIGH	200 IU/mL
REAGENT BLANK	REAG/DIL	NUMBER OF STEPS		1
CLEANER	NO	CALCULATION STEPS A		ENDPOINT
WAVELENGTH	550 nm	READING FIRST		2
DECIMAL POSITION	1	LAST		7
UNIT	IU/mL			
ANALYSIS		CALIBRATION		
POST DIL. FACTOR	2	CALIB. INTERVAL		ON REQUEST
CONC. FACTOR	NO	TIME		NO
SAMPLE CYCLE	1	BLANK RANGE	LOW	NO
VOLUME	3,0 µL		HIGH	NO
DILUTION NAME	H2O	REAGENT RANGE	LOW	NO
VOLUME	20,0 µL		HIGH	NO
REAGENT CYCLE	1	STANDARD POS		*
VOLUME	300 µL	STD-1		@
		STD-2		
		STD-3		
CALCULATION				
SAMPLE LIMIT	NO			
REACTION DIRECTION	INCREASE	REPLICATE		DUPL
CHECK	ON	DEVIATION		15%
CONVERSION FACTOR	1.000	CS1 POS		*
OFFSET	0.000	CS2 POS		*
		CS3 POS		*

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GENERAL		CONVERSION FACTOR	
MEASUREMENT MODE	ABSORB	OFFSET	1.000
REACTION MODE	D-R-S	TEST RANGE LOW	0.000
CALIBRATION MODE	SLOP AVG	HIGH	20 mg/dL
REAGENT BLANK	REAG/DIL	NORM RANGE LOW	500 mg/dL
CLEANER	NO	HIGH	70 mg/dL
WAVELENGTH	340 nm	NUMBER OF STEPS	150 mg/dL
DECIMAL POSITION	1	CALCULATION STEPS A	1
UNIT	mg/dL	READING FIRST	ENDPOINT
		LAST	1
ANALYSIS			25
DILUTION NAME	NaCl	CALIBRATION	
FACTOR	5	CALIB. INTERVAL	ON REQUEST
TIME	NO	TIME	NO
STD	MAIN INDIRECT	BLANK RANGE LOW	NO
MAIN STD	##	HIGH	NO
POST DIL. FACTOR	2	REAGENT RANGE LOW	NO
CONC. FACTOR	NO	HIGH	NO
SAMPLE CYCLE	1	STANDARD POS	*
VOLUME	5 µL	STD-1	@
DILUTION NAME	NaCl	STD-2	
VOLUME	00 µL	STD-3	
REAGENT CYCLE	1	REPLICATE	DUPL
VOLUME	250 µL	DEVIATION	15%
CALCULATION		CONTROL	
SAMPLE LIMIT	NO	CS1 POS	*
REACTION DIRECTION	INCREASE	CS2 POS	*
CHECK	ON	CS3 POS	*

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GENERAL		CONVERSION FACTOR	
MEASUREMENT MODE	ABSORB	OFFSET	1.000
REACTION MODE	D-R-S	TEST RANGE LOW	0.000
CALIBRATION MODE	SLOP AVG	HIGH	5 mg/dL
REAGENT BLANK	REAG/DIL	NORM RANGE LOW	100 mg/dL
CLEANER	NO	HIGH	10 mg/dL
WAVELENGTH	340 nm	NUMBER OF STEPS	40 mg/dL
DECIMAL POSITION	1	CALCULATION STEPS A	1
UNIT	mg/dL	READING FIRST	ENDPOINT
		LAST	1
ANALYSIS			25
DILUTION NAME	NaCl	CALIBRATION	
FACTOR	5	CALIB. INTERVAL	ON REQUEST
TIME	NO	TIME	NO
STD	MAIN INDIRECT	BLANK RANGE LOW	NO
MAIN STD	##	HIGH	NO
POST DIL. FACTOR	2	REAGENT RANGE LOW	NO
CONC. FACTOR	NO	HIGH	NO
SAMPLE CYCLE	1	STANDARD POS	*
VOLUME	12,5 µL	STD-1	@
DILUTION NAME	NaCl	STD-2	
VOLUME	00 µL	STD-3	
REAGENT CYCLE	1	REPLICATE	DUPL
VOLUME	250 µL	DEVIATION	15%
CALCULATION		CONTROL	
SAMPLE LIMIT	NO	CS1 POS	*
REACTION DIRECTION	INCREASE	CS2 POS	*
CHECK	ON	CS3 POS	*

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GENERAL		SAMPLE LIMIT	NO
MEASUREMENT MODE	ABSORB	REACTION DIRECTION	INCREASE
REACTION MODE	D-R-S-SR1	CHECK	ON
CALIBRATION MODE	LOGIT/LOG5	CONVERSION FACTOR	1.000
REAGENT BLANK	REAG/DIL	OFFSET	0.000
CLEANER	NO	TEST RANGE LOW	ON
WAVELENGTH	600 nm	HIGH	ON
DECIMAL POSITION	2	NORM RANGE LOW	20 µg/L
UNIT	µg/L	HIGH	250 µg/L
ANALYSIS		NUMBER OF STEPS	1
DILUTION NAME	H2O	CALCULATION STEPS A	ENDPOINT
FACTOR	NO	READING FIRST	3
TIME	NO	LAST	23
STD	MAIN INDIRECT	CALIBRATION	
MAIN STD	@	CALIB. INTERVAL	ON REQUEST
POST DIL. FACTOR	2	TIME	NO
CONC. FACTOR	NO	BLANK RANGE LOW	NO
SAMPLE CYCLE	2	REAGENT RANGE LOW	NO
VOLUME	28 µL	HIGH	NO
DILUTION NAME	H2O	STANDARD POS	*
VOLUME	5 µL	STD-1 - @	STD-2 - @÷2
REAGENT CYCLE	1	STD-3 - @÷4	STD-4 - @÷8
VOLUME	220 µL	STD-5 - @÷16	STD-6 - @÷32
START R1 CYCLE	1	REPLICATE	SINGLE
VOLUME	55	DEVIATION	NO
DILUTION NAME	H2O	CONTROL	
VOLUME	5 µL	CS1 POS	*

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GENERAL		SAMPLE LIMIT	NO
MEASUREMENT MODE	ABSORB	REACTION DIRECTION	INCREASE
REACTION MODE	D-R-S-SR1	CHECK	ON
CALIBRATION MODE	LOGIT/LOG5	CONVERSION FACTOR	1.000
REAGENT BLANK	REAG/DIL	OFFSET	0.000
CLEANER	NO	TEST RANGE LOW	ON
WAVELENGTH	600 nm	HIGH	ON
DECIMAL POSITION	2	NORM RANGE LOW	NO
UNIT	IU/mL	HIGH	30 IU/mL
ANALYSIS		NUMBER OF STEPS	1
DILUTION NAME	H2O	CALCULATION STEPS A	ENDPOINT
FACTOR	NO	READING FIRST	1
TIME	NO	LAST	7
STD	MAIN DIRECT	CALIBRATION	
MAIN STD	@	CALIB. INTERVAL	ON REQUEST
POS	*	TIME	NO
POST DIL. FACTOR	2	BLANK RANGE LOW	NO
CONC. FACTOR	NO	REAGENT RANGE LOW	NO
SAMPLE CYCLE	2	HIGH	NO
VOLUME	2 µL	STANDARD POS	*
DILUTION NAME	H2O	STD-1 - @	STD-2 - @÷2
VOLUME	10 µL	STD-3 - @÷4	STD-4 - @÷8
REAGENT CYCLE	1	STD-5 - @÷16	STD-6 - @÷32
VOLUME	220 µL	STD-7 - 0.00	STD-8 - NO
START R1 CYCLE	1	REPLICATE	SINGLE
VOLUME	60	DEVIATION	NO
DILUTION NAME	H2O	CONTROL	
VOLUME	10 µL	CS1 POS	*

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GENERAL		TEST RANGE LOW	ON
MEASUREMENT MODE	ABSORB	HIGH	ON
REACTION MODE	R-S-SR1	NORM RANGE LOW	NO
CALIBRATION MODE	LOGIT/LOG5	HIGH	NO
REAGENT BLANK	NO	NUMBER OF STEPS	1
CLEANER	NO	CALCULATION STEPS A	ENDPOINT
WAVELENGTH	600 nm	READING FIRST	13
DECIMAL POSITION	2	LAST	24
UNIT	%		
ANALYSIS		CALIBRATION	
POST DIL. FACTOR	NO	CALIB. INTERVAL	ON REQUEST
CONC. FACTOR	NO	TIME	NO
SAMPLE CYCLE	1	BLANK RANGE LOW	NO
VOLUME	4 µL	HIGH	NO
DILUTION NAME	H2O	REAGENT RANGE LOW	NO
VOLUME	8 µL	HIGH	NO
REAGENT CYCLE	1	STANDARD POS	*
VOLUME	150 µL	STD-1	@
REAGENT CYCLE	12	STD-2	@
VOLUME	50 µL	STD-3	@
DIL	8 µL	STD-4	@
CALCULATION		STD-5	@
SAMPLE LIMIT	NO	REPLICATE	DUPL
REACTION DIRECTION	INCREASE	DEVIATION	10.0 %
CHECK	ON	CALC MODEL	Logit/LOG
CONVERSION FACTOR	1.000	CONTROL	
OFFSET	0.000	CS1 POS	*
		CS2 POS	*

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GENERAL		CONVERSION FACTOR	
MEASUREMENT MODE	ABSORB	OFFSET	1.000
REACTION MODE	D-R-S	TEST RANGE LOW	0.000
CALIBRATION MODE	SLOP AVG	HIGH	30 mg/dL
REAGENT BLANK	REAG/DIL	NORM RANGE LOW	600 mg/dL
CLEANER	NO	HIGH	40 mg/dL
WAVELENGTH	340 nm	NUMBER OF STEPS	390 mg/dL
DECIMAL POSITION	1	CALCULATION STEPS A	1
UNIT	mg/dL	READING FIRST	1
ANALYSIS		LAST	25
DILUTION NAME	NaCl	CALIBRATION	
FACTOR	5	CALIB. INTERVAL	ON REQUEST
TIME	NO	TIME	NO
STD	MAIN INDIRECT	BLANK RANGE LOW	NO
MAIN STD	##	HIGH	NO
POST DIL. FACTOR	2	REAGENT RANGE LOW	NO
CONC. FACTOR	NO	HIGH	NO
SAMPLE CYCLE	1	STANDARD POS	*
VOLUME	2.5 µL	STD-1	@
DILUTION NAME	NaCl	STD-2	
VOLUME	00 µL	STD-3	
REAGENT CYCLE	1	REPLICATE	DUPL
VOLUME	250 µL	DEVIATION	15%
CALCULATION		CONTROL	
SAMPLE LIMIT	NO	CS1 POS	*
REACTION DIRECTION	INCREASE	CS2 POS	*
CHECK	ON	CS3 POS	*

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GENERAL		CONVERSION FACTOR	
MEASUREMENT MODE	ABSORB	OFFSET	1.000
REACTION MODE	D-R-S	TEST RANGE LOW	0.000
CALIBRATION MODE	SLOP AVG	HIGH	250 mg/dL
REAGENT BLANK	REAG/DIL	NORM RANGE LOW	3000 mg/dL
CLEANER	NO	HIGH	525 mg/dL
WAVELENGTH	340 nm	NUMBER OF STEPS	1650 mg/dL
DECIMAL POSITION	1	CALCULATION STEPS A	1
UNIT	mg/dL	READING FIRST	ENDPOINT
ANALYSIS		LAST	25
DILUTION NAME	NaCl	CALIBRATION	
FACTOR	15	CALIB. INTERVAL	ON REQUEST
TIME	NO	TIME	NO
STD	MAIN INDIRECT	BLANK RANGE LOW	NO
MAIN STD	##	HIGH	NO
POST DIL. FACTOR	2	REAGENT RANGE LOW	NO
CONC. FACTOR	NO	HIGH	NO
SAMPLE CYCLE	1	STANDARD POS	*
VOLUME	2.5 µL	STD-1	@
DILUTION NAME	NaCl	STD-2	
VOLUME	00 µL	STD-3	
REAGENT CYCLE	1	REPLICATE	DUPL
VOLUME	250 µL	DEVIATION	15%
CALCULATION		CONTROL	
SAMPLE LIMIT	NO	CS1 POS	*
REACTION DIRECTION	INCREASE	CS2 POS	*
CHECK	ON	CS3 POS	*

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GENERAL		CONVERSION FACTOR	
MEASUREMENT MODE	ABSORB	OFFSET	1.000
REACTION MODE	D-R-S	TEST RANGE LOW	0.000
CALIBRATION MODE	SLOP AVG	HIGH	0 mg/dL
REAGENT BLANK	REAG/DIL	NORM RANGE LOW	500 mg/dL
CLEANER	NO	HIGH	25 mg/dL
WAVELENGTH	340 nm	NUMBER OF STEPS	310 mg/dL
DECIMAL POSITION	1	CALCULATION STEPS A	1
UNIT	mg/dL	READING FIRST	ENDPOINT
ANALYSIS		LAST	25
DILUTION NAME	NaCl	CALIBRATION	
FACTOR	10	CALIB. INTERVAL	ON REQUEST
TIME	NO	TIME	NO
STD	MAIN INDIRECT	BLANK RANGE LOW	NO
MAIN STD	##	HIGH	NO
POST DIL. FACTOR	2	REAGENT RANGE LOW	NO
CONC. FACTOR	NO	HIGH	NO
SAMPLE CYCLE	1	STANDARD POS	*
VOLUME	12.5 µL	STD-1	@
DILUTION NAME	NaCl	STD-2	
VOLUME	00 µL	STD-3	
REAGENT CYCLE	1	REPLICATE	DUPL
VOLUME	250 µL	DEVIATION	15%
CALCULATION		CONTROL	
SAMPLE LIMIT	NO	CS1 POS	*
REACTION DIRECTION	INCREASE	CS2 POS	*
CHECK	ON	CS3 POS	*

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GENERAL		CALCULATION			
MEASUREMENT MODE	ABSORB	SAMPLE LIMIT		NO	
REACTION MODE	D-R-S-SR1	REACTION DIRECTION		INCREASE	
CALIBRATION MODE	LOGIT/LOG5	CHECK		ON	
REAGENT BLANK	REAG/DIL	CONVERSION FACTOR		1.000	
CLEANER	NO	OFFSET		0.000	
WAVELENGTH	550 nm	TEST RANGE LOW		ON	
DECIMAL POSITION	1	HIGH		ON	
UNIT	mg/dL	NORM RANGE LOW		0.0 mg/dL	
		HIGH		30.0 mg/dL	
ANALYSIS		NUMBERS OF STEPS		1	
DILUTION NAME	H2O	CALCULATION STEPS A		ENDPOINT	
FACTOR	NO	READING FIRST		2	
TIME	NO	LAST		20	
STD	MAIN INDIRECT				
MAIN STD	@	CALIBRATION			
POST DIL. FACTOR	2.00	CALIB. INTERVAL		ON REQUEST	
CONC. FACTOR	NO	BLANK RANGE	LOW	NO	
SAMPLE CYCLE	3	BLANK RANGE	HIGH	NO	
VOLUME	5 µL	LOW		NO	
DILUTION NAME	H2O	REAGENT RANGE HIGH		NO	
VOLUME	10.0 µL	STANDARD POS		*	
REAGENT CYCLE	1	STD1	@	STD2	@÷2
VOLUME	225 µL	STD3	@÷4	STD4	@÷8
START R1 CYCLE	1	STD5	@÷16	STD6	0.0
VOLUME	60 µL	REPLICATE		SINGLE	
DILUTION VOLUME	10 µL	DEVIATION		15%	
CALCULATION		CONTROL			
SAMPLE LIMIT	NO	CS1 POS			
REACTION DIRECTION	INCREASE	CS2 POS		*	
CHECK	ON	CS3 POS		*	

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GENERAL		TEST RANGE	LOW	2	mg/L
MEASUREMENT MODE	ABSORB		HIGH	160	mg/L
REACTION MODE	R-S-SR1	NORM RANGE	LOW	0	mg/L
CALIBRATION MODE	SLOP AVG		HIGH	30	mg/L
REAGENT BLANK	REAG/DIL	NUMBER OF STEPS		1	
CLEANER	NO	CALCULATION STEPS A		ENDPOINT	
WAVELENGTH	550 nm	READING FIRST		2	
DECIMAL POSITION	1	LAST		7	
UNIT	mg/L				
ANALYSIS		CALIBRATION			
POST DIL. FACTOR	2	CALIB. INTERVAL		ON REQUEST	
CONC. FACTOR	NO	TIME		NO	
SAMPLE CYCLE	1	BLANK RANGE	LOW	NO	
VOLUME	2.0 µL		HIGH	NO	
DILUTION NAME	H2O	REAGENT RANGE	LOW	NO	
VOLUME	10.0 µL		HIGH	NO	
REAGENT CYCLE	1	STANDARD POS		*	
VOLUME	200 µL	STD-1		@	
START R1 CYCLE	2	STD-2			
VOLUME	50 µL	STD-3			
DILUTION VOLUME	0	REPLICATE		DUPL	
		DEVIATION		15%	
CALCULATION		CONTROL			
SAMPLE LIMIT	NO	CS1 POS		*	
REACTION DIRECTION	INCREASE	CS2 POS		*	
CHECK	ON	CS3 POS		*	
CONVERSION FACTOR	1.000				
OFFSET	0.000				

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GENERAL		TEST RANGE	LOW	0 mg/L
MEASUREMENT MODE	ABSORB		HIGH	150 mg/L
REACTION MODE	R-S	NORM RANGE	LOW	NO
CALIBRATION MODE	SLOP AVG		HIGH	NO
REAGENT BLANK	REAG/DIL	NUMBER OF STEPS		1
CLEANER	NO	CALCULATION STEPS A		ENDPOINT
WAVELENGTH	550 nm	READING FIRST		2
DECIMAL POSITION	1	LAST		6
UNIT	mg/L			
ANALYSIS		CALIBRATION		
POST DIL. FACTOR	2	CALIB. INTERVAL		ON REQUEST
CONC. FACTOR	NO	TIME		NO
SAMPLE CYCLE	1	BLANK RANGE	LOW	NO
VOLUME	2 µL		HIGH	NO
DILUTION NAME	H2O	REAGENT RANGE	LOW	NO
VOLUME	10.0 µL		HIGH	NO
REAGENT CYCLE	1	STANDARD POS		*
VOLUME	300 µL	STD-1		@
		STD-2		
		STD-3		
CALCULATION		CONTROL		
SAMPLE LIMIT	NO	REPLICATE		DUPL
REACTION DIRECTION	INCREASE	DEVIATION		15%
CHECK	ON			
CONVERSION FACTOR	1.000	CS1 POS		*
OFFSET	0.000	CS2 POS		*
		CS3 POS		*

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GENERAL		OFFSET		0.000
MEASUREMENT MODE	ABSORB	TEST RANGE	LOW	0
REACTION MODE	R-S-SR1		HIGH	150
CALIBRATION MODE	SLOPE AVG	NORM RANGE	LOW	0 mg/L
REAGENT BLANK	REAG/DIL		HIGH	6 mg/L
CLEANER	NO	NUMBER OF STEPS		1
WAVELENGTH	550 nm	CALCULATION STEPS A		ENDPOINT
DECIMAL POSITION	1	READING FIRST		2
UNIT	mg/L		LAST	7
ANALYSIS		CALIBRATION		
POS DIL. FACTOR	2	CALIB. INTERVAL		ON REQUEST
CONC. FACTOR	NO	REAGENT RANGE	LOW	NO
SAMPLE CYCLE	2		HIGH	NO
VOLUME	2,0 µL	BLANK RANGE	LOW	NO
DILUTION NAME	H2O		HIGH	NO
VOLUME	10.0 µL	STANDARD POS		*
REAGENT CYCLE	1	STD 1:		@
VOLUME	300 µL	STD 2:		
START R1	1	STD 3:		
VOLUME	75 µL			
DILUENT NAME	H2O	REPLICATE		DUPL
VOLUME	0.0	DEVIATION		15%
		CORRECTION STD		NO
CALCULATION		CONTROL		
SAMPLE LIMIT	NO	CS1 POS		*
REACTION DIRECTION	INCREASE	CS2 POS		*
CHECK	ON	CS3 POS		*
CONVERS FACTOR	1.000			

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GENERAL		OFFSET		0.000
MEASUREMENT MODE	ABSORB	TEST RANGE	LOW	NO
REACTION MODE	D-R-S		HIGH	NO
CALIBRATION MODE	LOGIT/LOG5	NORM RANGE	LOW	0 mg/L
REAGENT BLANK	REAG/DIL		HIGH	3 mg/L
CLEANER	NO	NUMBER OF STEPS		1
WAVELENGTH	550 nm	CALCULATION STEPS A		ENDPOINT
DECIMAL POSITION	2	READING FIRST		2
UNIT	mg/L		LAST	14
ANALYSIS		CALIBRATION		
DILUTION NAME	H2O	CALIB. INTERVAL		ON REQUEST
FACTOR	NO	REAGENT RANGE	LOW	NO
TIME	NO		HIGH	NO
STD	MAIN DIRECT	BLANK RANGE	LOW	NO
MAIN STD	@		HIGH	NO
POS DIL. FACTOR	2	STANDARD		POS *
CONC. FACTOR	NO	1: @		2: @÷2
SAMPLE CYCLE	2	3: @÷4		4: @÷8
VOLUME	2,5 µL	5: @÷16		6: 0.00
DILUTION NAME	H2O	7: NO		8: NO
VOLUME	10.0 µL	REPLICATE		DUPL
REAGENT CYCLE	1	DEVIATION		15%
VOLUME	250 µL	CORRECTION STD		NO
CALCULATION		CONTROL		
SAMPLE LIMIT	NO	CS1 POS		*
REACTION DIRECTION	INCREASE	CS2 POS		*
CHECK	ON	CS3 POS		*
CONVERS FACTOR	1.000			

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GENERAL		CONVERSION FACTOR	
MEASUREMENT MODE	ABSORB	OFFSET	1.000
REACTION MODE	D-R-S	TEST RANGE LOW	0.000
CALIBRATION MODE	SLOP AVG	HIGH	30 mg/dL
REAGENT BLANK	REAG/DIL	NORM RANGE LOW	500 mg/dL
CLEANER	NO	HIGH	200 mg/dL
WAVELENGTH	340 nm	NUMBER OF STEPS	350 mg/dL
DECIMAL POSITION	1	CALCULATION STEPS A	1
UNIT	mg/dL	READING FIRST	ENDPOINT
		LAST	1
ANALYSIS			25
DILUTION NAME	NaCl	CALIBRATION	
FACTOR	15	CALIB. INTERVAL	ON REQUEST
TIME	NO	TIME	NO
STD	MAIN INDIRECT	BLANK RANGE LOW	NO
MAIN STD	##	HIGH	NO
POST DIL. FACTOR	2	REAGENT RANGE LOW	NO
CONC. FACTOR	NO	HIGH	NO
SAMPLE CYCLE	1	STANDARD POS	*
VOLUME	5 µL	STD-1	@
DILUTION NAME	NaCl	STD-2	
VOLUME	00 µL	STD-3	
REAGENT CYCLE	1	REPLICATE	DUPL
VOLUME	250 µL	DEVIATION	15%
CALCULATION		CONTROL	
SAMPLE LIMIT	NO	CS1 POS	*
REACTION DIRECTION	INCREASE	CS2 POS	*
CHECK	ON	CS3 POS	*

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