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EXPRESS

550

BIOQUÍMICA

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Test Name:	ÁCIDO ÚRICO	Test:	ACU
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	510	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.10	High Blank A Limit:	0.30
Low A Limit:	0.00	High A Limit:	1.10
Low Normal:	*	High Normal:	*
Linearity Limit:	20	Curve S.D. Limit:	*
Test Name:	ÁCIDO ÚRICO	TEST:	ACU
Test Bar Code:			
Sample Volume:	5 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	250	ACU	490
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	ALBUMINA	Test:	ALB
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	g/dL	No. of Decimal Places:	2
Primary Wavelength:	600	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	0.50
Low A Limit:	-0.25	High A Limit:	0.90
Low Normal:	*	High Normal:	*
Linearity Limit:	6	Curve S.D. Limit:	*
Test Name:	ALBUMINA	TEST:	ALB
Test Bar Code:			
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	400	ALB	120
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	AMILASE	Test:	AMI
Test Bar Code:			
Test Type:	Kinetic	Curve Type:	Blanked Linear
Units:	U/L	No. of Decimal Places:	0
Primary Wavelength:	405	Secondary Wavelength	
Read Time/Interval:	60	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	1.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	1038	Curve S.D. Limit:	*
Test Name:	AMILASE	TEST:	AMI
Test Bar Code:			
Sample Volume:	3 μ L	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	AMI	60
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	ALT/TGP	Test:	ALT
Test Bar Code:			
Test Type:	Kinetic	Curve Type:	Blanked Linear
Units:	U/L	No. of Decimal Places:	0
Primary Wavelength:	340	Secondary Wavelength	
Read Time/Interval:	60	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.80	High Blank A Limit:	2.000
Low A Limit:	0.45	High A Limit:	2.000
Low Normal:	*	High Normal:	*
Linearity Limit:	350	Curve S.D. Limit:	*
Test Name:	ALT/TGP	TEST:	ALT
Test Bar Code:			
Sample Volume:	30 µL	Sample Diluent:	
Rerun Dilution Ratio:	5	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	ALT	60
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	AST/TGO	Test:	AST
Test Bar Code:			
Test Type:	Kinetic	Curve Type:	Blanked Linear
Units:	U/L	No. of Decimal Places:	0
Primary Wavelength:	340	Secondary Wavelength	
Read Time/Interval:	60	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.80	High Blank A Limit:	2.000
Low A Limit:	0.45	High A Limit:	2.000
Low Normal:	*	High Normal:	*
Linearity Limit:	440	Curve S.D. Limit:	*
Test Name:	AST/TGO	TEST:	AST
Test Bar Code:			
Sample Volume:	30 µL	Sample Diluent:	
Rerun Dilution Ratio:	5	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	AST	60
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	BIL. DIRETA	Test:	BD
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	2
Primary Wavelength:	540	Secondary Wavelength	None
Read Time/Interval:	10	Sample Blank:	Si
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.10	High Blank A Limit:	0.30
Low A Limit:	-0.10	High A Limit:	1.10
Low Normal:	*	High Normal:	*
Linearity Limit:	15.00	Curve S.D. Limit:	*
Test Name:	BIL. DIRETA	TEST:	BD
Test Bar Code:			
Sample Volume:	30 µL	Sample Diluent:	
Rerun Dilution Ratio:	5	Predilution Ratio:	1
		Reagent/Diluent:	BD
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	270	BD	30
REAGENT 2	20	NIT	300
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	CÁLCIO CPC	Test:	Ca
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	570	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.10	High Blank A Limit:	0.20
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	20	Curve S.D. Limit:	*
Test Name:	CÁLCIO CPC	TEST:	Ca
Test Bar Code:		Reagent Diluent:	CaA
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	150	CaB	150
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	CK NAC	Test:	CK
Test Bar Code:			
Test Type:	Kinetic	Curve Type:	Blanked Linear
Units:	U/L	No. of Decimal Places:	0
Primary Wavelength:	340	Secondary Wavelength	
Read Time/Interval:	60	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	2000	Curve S.D. Limit:	*
Test Name:	CK NAC	TEST:	CK
Test Bar Code:			
Sample Volume:	10 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	250	CK	120
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	CK MB	Test:	MB
Test Bar Code:			
Test Type:	Kinetic	Curve Type:	Blanked Linear
Units:	U/L	No. of Decimal Places:	0
Primary Wavelength:	340	Secondary Wavelength	
Read Time/Interval:	60	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	600	Curve S.D. Limit:	
Test Name:	CK MB	TEST:	MB
Test Bar Code:			
Sample Volume:	10 µL	Sample Diluent:	
Rerun Dilution Ratio:	2	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	250	MB	LAG TIME
REAGENT 2			180
REAGENT 3			
REAGENT 4			
CONTROLS:			

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TEST NAME:	CLORETOS	TEST:	CL
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mEq/L	No. of Decimal Places:	0
Primary Wavelength:	510	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	0.50
Low A Limit:	0.00	High A Limit:	0.50
Low Normal:	*	High Normal:	*
Linearity Limit:	150	Curve S.D. Limit:	*
Test Name:	CLORETOS	TEST:	CL
Test Bar Code:			
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	400	CL	290
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	COLESTEROL	Test:	COL
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	0
Primary Wavelength:	510	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.10	High Blank A Limit:	0.30
Low A Limit:	-0.10	High A Limit:	1.00
Low Normal:	*	High Normal:	*
Linearity Limit:	800	Curve S.D. Limit:	*
Test Name:	COLESTEROL	TEST:	COL
Test Bar Code:			
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	COL	LAG TIME
REAGENT 2			300
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	CREATININA	Test:	CRE
Test Bar Code:			
Test Type:	Two Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	2
Primary Wavelength:	510	Secondary Wavelength	
Read Time/Interval:	60	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	0
Low Blank A Limit:	0.00	High Blank A Limit:	1.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	12	Curve S.D. Limit:	*
Test Name:	CREATININA	TEST:	CRE
Test Bar Code:			
Sample Volume:	25 µL	Sample Diluent:	
Rerun Dilution Ratio:	5	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	250	CRE	LAG TIME
REAGENT 2			30
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	DESIDROGENASE	Test:	LDH
Test Bar Code:			
Test Type:	Kinetic	Curve Type:	Blanked Linear
Units:	U/L	No. of Decimal Places:	0
Primary Wavelength:	340	Secondary Wavelength	
Read Time/Interval:	60	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. of Replicates:	2
Low Blank A Limit:	0.80	High Blank A Limit:	2.00
Low A Limit:	0.45	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	2000	Curve S.D. Limit:	*
Test Name:	DESIDROGENASE	TEST:	LDH
Test Bar Code:			
Sample Volume:	5 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	250	LDH	120
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	FERRO FEZ	Test:	Fe
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	µg/dL	No. of Decimal Places:	0
Primary Wavelength:	570	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	Yes
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	0.80
Low A Limit:	0.01	High A Limit:	0.80
Low Normal:	*	High Normal:	*
Linearity Limit:	500	Curve S.D. Limit:	*
Test Name:	FERRO FEZ	TEST:	Fe
Test Bar Code:			
Sample Volume:	30 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	250	Fe1	290
REAGENT 2	60	Fe2	300
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	FOSF. ALCALINA	Test:	FAL
Test Bar Code:			
Test Type:	Kinetic	Curve Type:	Blanked Linear
Units:	U/L	No. of Decimal Places:	0
Primary Wavelength:	405	Secondary Wavelength	
Read Time/Interval:	60	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.01	High Blank A Limit:	1.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	700	Curve S.D. Limit:	*
Test Name:	FOSF. ALCALINA	TEST:	FAL
Test Bar Code:			
Sample Volume:	5 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	250	FAL	60
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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TEST NAME:	FÓSFORO UV	TEST:	FUV
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	340	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	No
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.10	High Blank A Limit:	1.50
Low A Limit:	0.00	High A Limit:	1.50
Low Normal:	*	High Normal:	*
Linearity Limit:	15	Curve S.D. Limit:	*
Test Name:	FÓSFORO UV	TEST:	FUV
Test Bar Code:			
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	300	FUV	300
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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TEST NAME:	FRUTOSAMINA	TEST:	FRU
Test Bar Code:			
Test Type:	Two Point	Curve Type:	Blanked Linear
Units:	mmol/L	No. of Decimal Places:	1
Primary Wavelength:	546	Secondary Wavelength	
Read Time/Interval:	180	Sample Blank:	No
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.10	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	8.5	Curve S.D. Limit:	*
Test Name:	FRUTOSAMINA	TEST:	FRU
Test Bar Code:			
Sample Volume:	12 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	240	FRU	400
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	GAMA GT	Test:	GGT
Test Bar Code:			
Test Type:	Kinetic	Curve Type:	Blanked Linear
Units:	U/L	No. of Decimal Places:	0
Primary Wavelength:	405	Secondary Wavelength	
Read Time/Interval:	60	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.01	High Blank A Limit:	1.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	250	Curve S.D. Limit:	*
Test Name:	GAMA GT	TEST:	GGT
Test Bar Code:			
Sample Volume:	12 µL	Sample Diluent:	
Rerun Dilution Ratio:	2	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	250	GGT	
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			
			LAG TIME
			60

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Test Name:	GLICOSE	Test:	GLI
Test Bar Code:			
Test Type:	Two Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	0
Primary Wavelength:	510	Secondary Wavelength	
Read Time/Interval:	60	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.30	High Blank A Limit:	0.30
Low A Limit:	-0.10	High A Limit:	1.50
Low Normal:	*	High Normal:	*
Linearity Limit:	400	Curve S.D. Limit:	*
Test Name:	GLICOSE	TEST:	GLI
Test Bar Code:			
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	300	GLI	30
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	HDL COLESTEROL	Test:	HDL
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	0
Primary Wavelength:	510	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.10	High Blank A Limit:	0.30
Low A Limit:	-0.10	High A Limit:	1.00
Low Normal:	*	High Normal:	*
Linearity Limit:	150	Curve S.D. Limit:	*
Test Name:	HDL COLESTEROL	TEST:	HDL COL
Test Bar Code:			
Sample Volume:	25 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	250	HDL COL	490
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	HDL	Test:	HDL
Test Bar Code:			
Test Type:	Two Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	0
Primary Wavelength:	600	Secondary Wavelength	
Read Time/Interval:	290	Sample Blank:	No
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.01	High Blank A Limit:	0.50
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	180	Curve S.D. Limit:	*
Test Name:	HDL	TEST:	HDL
Test Bar Code:			
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	188	HDL1	300
REAGENT 2	62	HDL2	10
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	HDL	Test:	HDL
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	0
Primary Wavelength:	600	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	Sim
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.01	High Blank A Limit:	0.50
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	180	Curve S.D. Limit:	*
Test Name:	HDL	TEST:	HDL
Test Bar Code:			
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	250	HDL1	300
REAGENT 2	80	HDL2	270
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	LACTATO	Test:	LACT
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	546	Secondary Wavelength	
Read Time/Interval:	5	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.10	High Blank A Limit:	0.30
Low A Limit:	0.00	High A Limit:	1.10
Low Normal:	*	High Normal:	*
Linearity Limit:	120	Curve S.D. Limit:	*
Test Name:	LACTATO	TEST:	LACT
Test Bar Code:			
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	LACT	300
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	LDL	Test:	LDL
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	0
Primary Wavelength:	600	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	Sim
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.01	High Blank A Limit:	0.50
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	400	Curve S.D. Limit:	*
Test Name:	LDL	TEST:	LDL
Test Bar Code:	LDL		
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	250	LDL1	300
REAGENT 2	80	LDL2	270
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	MAGNÉSIO	Test:	MG
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	510	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	*	High Blank A Limit:	*
Low A Limit:	*	High A Limit:	*
Low Normal:	*	High Normal:	*
Linearity Limit:	4.5	Curve S.D. Limit:	*
Test Name:	MAGNÉSIO	TEST:	MG
Test Bar Code:			
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	MG	LAG TIME
REAGENT 2			120
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	PROTEÍNA TOTAL	Test:	PRO
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	g/dL	No. of Decimal Places:	2
Primary Wavelength:	540	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.01	High Blank A Limit:	1.00
Low A Limit:	-0.25	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	12	Curve S.D. Limit:	*
Test Name:	PROTEÍNA	TEST:	PRO
Test Bar Code:	PRO		
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	PRO	590
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	PROTEÍNA URINÁRIA	Test:	PRO-U
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/L	No. of Decimal Places:	0
Primary Wavelength:	600	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	0.50
Low A Limit:	0.00	High A Limit:	1.50
Low Normal:	*	High Normal:	*
Linearity Limit:	4000	Curve S.D. Limit:	*
Test Name:	PROTEÍNA	TEST:	PRO-U
Test Bar Code:			
Sample Volume:	6 µL	Sample Diluent:	
Rerun Dilution Ratio:	2	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	300	PRO-U	290
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	TRIGLICÉRIDES	Test:	TGR
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	0
Primary Wavelength:	510	Secondary Wavelength	
Read Time/Interval:	10	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	-0.01	High Blank A Limit:	0.30
Low A Limit:	-0.01	High A Limit:	1.60
Low Normal:	*	High Normal:	*
Linearity Limit:	1000	Curve S.D. Limit:	*
Test Name:	TRIGLICERIDES	TEST:	TGR
Test Bar Code:			
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	300	TGR	490
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	URÉIA UV	Test:	UUV
Test Bar Code:			
Test Type:	Kinetic	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	0
Primary Wavelength:	340	Secondary Wavelength	
Read Time/Interval:	60	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	1.00	High Blank A Limit:	2.00
Low A Limit:	0.60	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	250	Curve S.D. Limit:	*
Test Name:	URÉIA UV	TEST:	UUV
Test Bar Code:			
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	UUV	LAG TIME
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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