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EXPRESS

550

TURBIDIMETRIA

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Test Name:	ALFAGLICO	Test:	AGP
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	340	Secondary Wavelength	None
Read Time/Interval:	20	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	200	Curve S.D. Limit:	*
Test Name:	ALFAGLICO	TEST:	AGP
Test Bar Code:			
Sample Volume:	15 μ L	Sample Diluent:	NaCl
Rerun Dilution Ratio:	1	Predilution Ratio:	10
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	300	AGP	360
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	APO-A	Test:	APO-A
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	340	Secondary Wavelength	None
Read Time/Interval:	20	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	300	Curve S.D. Limit:	*
Test Name:	APO-A	TEST:	APO-A
Test Bar Code:			
Sample Volume:	6 µL	Sample Diluent:	NaCl
Rerun Dilution Ratio:	2	Predilution Ratio:	10
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	APOA	360
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	APO-B	Test:	APO-B
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	340	Secondary Wavelength	None
Read Time/Interval:	20	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	300	Curve S.D. Limit:	*
Test Name:	APO-B	TEST:	APO-B
Test Bar Code:			
Sample Volume:	7 µL	Sample Diluent:	NaCl
Rerun Dilution Ratio:	2	Predilution Ratio:	10
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	280	APO-B	360
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	ASO	Test:	ASO
Test Bar Code:			
Test Type:	Two Point	Curve Type:	Blanked Linear
Units:	UI/mL	No. of Decimal Places:	1
Primary Wavelength:	540	Secondary Wavelength	None
Read Time/Interval:	120	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	800	Curve S.D. Limit:	*
Test Name:	ASO	TEST:	ASOD
Test Bar Code:	ASO		
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	30	ASOL	270
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	COMPLEMENTO C3	Test:	C3
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	340	Secondary Wavelength	None
Read Time/Interval:	20	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.000	High Blank A Limit:	2.00
Low A Limit:	0.000	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	500	Curve S.D. Limit:	*
Test Name:	COMPLEMENTO C3	TEST:	C3
Test Bar Code:			
Sample Volume:	6 µL	Sample Diluent:	NaCl
Rerun Dilution Ratio:	1	Predilution Ratio:	10
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	C3	360
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	COMPLEMENTO C4	Test:	C4
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	340	Secondary Wavelength	None
Read Time/Interval:	20	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	100	Curve S.D. Limit:	*
Test Name:	COMPLEMENTO C4	TEST:	C4
Test Bar Code:			
Sample Volume:	15 µL	Sample Diluent:	NaCl
Rerun Dilution Ratio:	2	Predilution Ratio:	10
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	300	C4	360
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	FATOR REUMATÓIDE	Test:	FR
Test Bar Code:			
Test Type:	Two Point	Curve Type:	2-LOGIT5
Units:	IU/mL	No. of Decimal Places:	1
Primary Wavelength:	600	Secondary Wavelength	None
Read Time/Interval:	120	Sample Blank:	NO
Factor:	1.000		
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	5	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	100	Curve S.D. Limit:	*
Test Name:	FR	TEST:	FR
Test Bar Code:		Reagent Diluent:	FRD
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1		
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	60	FRL	240
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	HBA1C	Test:	HBA1C
Test Bar Code:			
Test Type:	Two Point	Curve Type:	2-LOGIT5
Units:	%	No. of Decimal Places:	1
Primary Wavelength:	600	Secondary Wavelength	None
Read Time/Interval:	290	Sample Blank:	NO
Factor:	1.000		
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	5	No. Of Replicates:	1
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:		Curve S.D. Limit:	*
Test Name:	HBA1C	TEST:	HbA1c
Test Bar Code:			
Sample Volume:	5 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	188	HBA1c1	300
REAGENT 2	62	HbA1c2	10
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	IgA	Test:	IgA
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	340	Secondary Wavelength	None
Read Time/Interval:	20	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	600	Curve S.D. Limit:	*
Test Name:	IgA	TEST:	IgA
Test Bar Code:			
Sample Volume:	6 µL	Sample Diluent:	NaCl
Rerun Dilution Ratio:	1	Predilution Ratio:	10
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	300	IgA	360
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	IgG	Test:	IgG
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	0
Primary Wavelength:	340	Secondary Wavelength	None
Read Time/Interval:	20	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	3000	Curve S.D. Limit:	*
Test Name:	IgG	TEST:	IgG
Test Bar Code:			
Sample Volume:	4 µL	Sample Diluent:	NaCl
Rerun Dilution Ratio:	1	Predilution Ratio:	20
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	300	IgG	360
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	IgM	Test:	IgM
Test Bar Code:			
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	340	Secondary Wavelength	None
Read Time/Interval:	20	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	500	Curve S.D. Limit:	*
Test Name:	IgM	TEST:	IgM
Test Bar Code:			
Sample Volume:	15 µL	Sample Diluent:	NaCl
Rerun Dilution Ratio:	1	Predilution Ratio:	10
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	IgM	360
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	Lpa	Test:	Lpa
Test Bar Code:			
Test Type:	Two Point	Curve Type:	2-LOGIT5
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	540	Secondary Wavelength	None
Read Time/Interval:	240	Sample Blank:	NO
Factor:	1.00		
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	5	No. Of Replicates:	1
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	80	Curve S.D. Limit:	*
Test Name:	Lpa	TEST:	Lpa
Test Bar Code:		Reagent Diluent:	LpaD
Sample Volume:	4 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	50	LpaL	200
REAGENT 2			10
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	MICROALBUMINÚRIA	Test:	MAL
Test Bar Code:			
Test Type:	Two Point	Curve Type:	Blanked Linear
Units:	mg/L	No. of Decimal Places:	2
Primary Wavelength:	540	Secondary Wavelength	None
Read Time/Interval:	120	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	80	Curve S.D. Limit:	*
Test Name:	MAL	TEST:	MAL
Test Bar Code:		Reagent Diluent:	MALD
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	1	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
			LAG TIME
REAGENT 1	60	MALL	240
REAGENT 2			10
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	PCR	Test:	PCR
Test Bar Code:			
Test Type:	Two Point	Curve Type:	Blanked Linear
Units:	mg/L	No. of Decimal Places:	1
Primary Wavelength:	540	Secondary Wavelength	None
Read Time/Interval:	120	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.00
Low A Limit:	0.00	High A Limit:	2.00
Low Normal:	*	High Normal:	*
Linearity Limit:	150	Curve S.D. Limit:	*
Test Name:	PCR	TEST:	PCR
Test Bar Code:			
Sample Volume:	3 µL	Sample Diluent:	
Rerun Dilution Ratio:	2	Predilution Ratio:	1
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	PCR	10
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	PCR-U	Test:	PCR-U
Test Bar Code:	PCR-U		
Test Type:	Two Point	Curve Type:	2-LOGIT5
Units:	mg/L	No. of Decimal Places:	1
Primary Wavelength:	540	Secondary Wavelength	None
Read Time/Interval:	240	Sample Blank:	NO
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	6	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.50
Low A Limit:	0.00	High A Limit:	2.50
Low Normal:	*	High Normal:	*
Linearity Limit:	5	Curve S.D. Limit:	*
Test Name:	PCR-U	TEST:	PCR-U
Test Bar Code:	PCR-U	Sample Diluent:	
Sample Volume:	3 µL	Predilution Ratio:	1
Rerun Dilution Ratio:	1	Diluent Reagent	PCR-U-D
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	20	PCR-UL	280
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			

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Test Name:	TRANSFERRINA	Test:	TFR
Test Bar Code:	TFR		
Test Type:	End Point	Curve Type:	Blanked Linear
Units:	mg/dL	No. of Decimal Places:	1
Primary Wavelength:	340	Secondary Wavelength	None
Read Time/Interval:	20	Sample Blank:	NO
Factor:			
Calibration Interval:	*		
Normalization Interval:	*		
No. Of Calibrators:	2	No. Of Replicates:	2
Low Blank A Limit:	0.00	High Blank A Limit:	2.50
Low A Limit:	0.00	High A Limit:	2.50
Low Normal:	*	High Normal:	*
Linearity Limit:	500	Curve S.D. Limit:	*
Test Name:	TFR	TEST:	TFR
Test Bar Code:			
Sample Volume:	4 µL	Sample Diluent:	NaCl
Rerun Dilution Ratio:	1	Predilution Ratio:	10
REAGENT DILUENTE			
	REAGENT VOLUME	BAR CODE	DILUENT VOLUME
REAGENT 1	300	TFR	360
REAGENT 2			
REAGENT 3			
REAGENT 4			
CONTROLS:			