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## SELECTRA I - II FLEXOR / SPINLAB

BIOQUÍMICA

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NAME	ALBUMINA	PROZONE CHECK	No
ABBR. NAME	ALB	REF. MALE LOW	*
MODE	Endpoint	REF. MALE HIGH	*
WAVELENGTH	620 nm	REF. FEMALE LOW	*
UNITS	g/dL	REF. FEMALE HIGH	*
DECIMALS	2	REF. PED. LOW	*
LOW CONC.	0 g/dL	REF. PED. HIGH	*
HIGH CONC.	6 g/dL	REF. PANIC LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PANIC HIGH	*
REPEAT	1	CONTROL 1	*
NUMBER	1	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	ALBUMINA		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	300 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	301 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	3 µL		
RERUN VOLUME	2 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	0 µL	Cal. High Limit	*
RERUN VOLUME	0 µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	4.5 min		
<b>MONO MODE</b>			
NAME	ALBUMINA		
SAMPLE BLANK	No	LOW ABSORBANCE	0.000 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	300 µL	R. ABS. L LIMIT	1.100 Abs
RERUN VOLUME	248 µL	R. ABS. H. LIMIT	2.000 Abs
SAMPLE			
NORMAL VOLUME	3 µL	REAGENT BLANK	Yes (0.000)
RERUN VOLUME	3.0 µL	Cal. Low Limit	*
INCUBATION TIME	11.5 min	FACTOR	#

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NAME	ALT/TGP	RERUN VOLUME	4 µL
ABBR. NAME	ALT	DELAY, MIN. TIME	77, 230 sec
MODE	Kinetic	PROZONE CHECK	No
WAVELENGTH	340 nm	REF. MALE LOW	*
UNITS	U/L	REF. MALE HIGH	*
DECIMALS	0	REF. FEMALE LOW	*
LOW CONC.	0 U/L	REF. FEMALE HIGH	*
HIGH CONC.	350 U/L	REF. PED. LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PED. HIGH	*
REPEAT	1	REF. PANIC LOW	*
NUMBER	1	REF. PANIC HIGH	*
CONCENTRATION	*	CONTROL 1	*
INTERVAL	0 days	CONTROL 2	*
CUT OFF	No	CONTROL 3	*
<b>DUAL MODE</b>			
NAME	ALT/TGP	CORRELAT. FACTOR	1.000
SAMPLE BLANK	No	CORRELAT. OFFSET	0.00
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.600 Abs
NORMAL VOLUME	220 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	236 µL		
SAMPLE		R. ABS. L LIMIT	- 0.100 Abs
NORMAL VOLUME	20 µL	R. ABS. H. LIMIT	0.500 Abs
RERUN VOLUME	4 µL		
R2 BOTTLE	5 mL	REAGENT BLANK	Yes (0.001)
NORMAL VOLUME	0 µL	Cal. Low Limit	*
RERUN VOLUME	0 µL	Cal. High Limit	*
PREDILUTION	No	FACTOR	#
INCUBATION TIME	No		
<b>MONO MODE</b>			
NAME	ALT/TGP	LOW ABSORBANCE	0.600 Abs
SAMPLE BLANK	No	HIGH ABSORBANCE	3.000 Abs
R1 BOTTLE	25 mL	R. ABS. L LIMIT	1.100 Abs
NORMAL VOLUME	220 µL	R. ABS. H. LIMIT	2.000 Abs
RERUN VOLUME	236 µL		
SAMPLE		REAGENT BLANK	Yes (0.000)
NORMAL VOLUME	20 µL	Cal. Low Limit	*
		FACTOR	#

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NAME	BDI	RERUN VOLUME	2 µL
ABBR. NAME	BDI	INCUBATION TIME	11.5
MODE	End Point	PROZONE CHECK	No
WAVELENGTH	546 nm	REF. MALE LOW	*
UNITS	mg/dL	REF. MALE HIGH	*
DECIMALS	2	REF. FEMALE LOW	*
LOW CONC.	0 mg/dL	REF. FEMALE HIGH	*
HIGH CONC.	15 mg/dL	REF. PED. LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PED. HIGH	*
REPEAT	1	REF. PANIC LOW	*
NUMBER	1	REF. PANIC HIGH	*
CONCENTRATION	*	CONTROL 1	*
INTERVAL	*	CONTROL 2	*
CUT OFF	No	CONTROL 3	*
<b>DUAL MODE</b>		CORRELAT. FACTOR	1.000
NAME	BDI	CORRELAT. OFFSET	0.00
SAMPLE BLANK	yes	LOW ABSORBANCE	0.600 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	240 µL	R. ABS. L LIMIT	- 0.100 Abs
RERUN VOLUME	248 µL	R. ABS. H. LIMIT	0.500 Abs
SAMPLE		REAGENT BLANK	Yes (0.001)
NORMAL VOLUME	10 µL		
RERUN VOLUME	2 µL	Cal. Low Limit	*
R2 BOTTLE	5 mL	Cal. High Limit	*
NORMAL VOLUME	0 µL	FACTOR	#
RERUN VOLUME	0 µL		
PREDILUTION	No	LOW ABSORBANCE	0.600 Abs
INCUBATION TIME	4.5	HIGH ABSORBANCE	3.000 Abs
<b>MONO MODE</b>		R. ABS. L LIMIT	1.100 Abs
NAME	BDI	R. ABS. H. LIMIT	2.000 Abs
SAMPLE BLANK	Yes		
R1 BOTTLE	25 mL	REAGENT BLANK	Yes (0.000)
NORMAL VOLUME	240 µL	Cal. Low Limit	*
RERUN VOLUME	248 µL	FACTOR	#
SAMPLE			
NORMAL VOLUME	10 µL		

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NAME	BTO	RERUN VOLUME	2 µL
ABBR. NAME	BTO	INCUBATION TIME	11.5
MODE	End Point	PROZONE CHECK	No
WAVELENGTH	546 nm	REF. MALE LOW	*
UNITS	mg/dL	REF. MALE HIGH	*
DECIMALS	2	REF. FEMALE LOW	*
LOW CONC.	0 mg/dL	REF. FEMALE HIGH	*
HIGH CONC.	15 mg/dL	REF. PED. LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PED. HIGH	*
REPEAT	1	REF. PANIC LOW	*
NUMBER	1	REF. PANIC HIGH	*
CONCENTRATION	*	CONTROL 1	*
INTERVAL	*	CONTROL 2	*
CUT OFF	No	CONTROL 3	*
<b>DUAL MODE</b>		CORRELAT. FACTOR	1.000
NAME	BTO	CORRELAT. OFFSET	0.00
SAMPLE BLANK	Yes	LOW ABSORBANCE	0.600 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	240 µL	R. ABS. L LIMIT	- 0.100 Abs
RERUN VOLUME	248 µL	R. ABS. H. LIMIT	0.500 Abs
SAMPLE		REAGENT BLANK	Yes (0.001)
NORMAL VOLUME	10 µL	Cal. Low Limit	*
RERUN VOLUME	2 µL	Cal. High Limit	*
R2 BOTTLE	5 mL	FACTOR	#
NORMAL VOLUME	0 µL	LOW ABSORBANCE	0.600 Abs
RERUN VOLUME	0 µL	HIGH ABSORBANCE	3.000 Abs
PREDILUTION	No	R. ABS. L LIMIT	1.100 Abs
INCUBATION TIME	4.5	R. ABS. H. LIMIT	2.000 Abs
<b>MONO MODE</b>		REAGENT BLANK	Yes (0.000)
NAME	BTO	Cal. Low Limit	*
SAMPLE BLANK	Yes	FACTOR	#
R1 BOTTLE	25 mL		
NORMAL VOLUME	240 µL		
RERUN VOLUME	248 µL		
SAMPLE			
NORMAL VOLUME	10 µL		

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NAME	CÁLCIO ARSENAZO	PROZONE CHECK	No
ABBR. NAME	CAAZ	REF. MALE LOW	*
MODE	Endpoint	REF. MALE HIGH	*
WAVELENGTH	620 nm	REF. FEMALE LOW	*
UNITS	Mg/dL	REF. FEMALE HIGH	*
DECIMALS	2	REF. PED. LOW	*
LOW CONC.	0.00 mg/dL	REF. PED. HIGH	*
HIGH CONC.	20 mg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PANIC HIGH	*
REPEAT	1	CONTROL 1	*
NUMBER	1	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	CÁLCIO ARSENAZO		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	240 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	242 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	4 µL		
RERUN VOLUME	2 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	0µL	Cal. High Limit	*
RERUN VOLUME	0µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	4.5 min		
<b>MONO MODE</b>			
NAME	CÁLCIO ARSENAZO		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	240 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	242 µL	R. ABS. L LIMIT	1.100 Abs
SAMPLE		R. ABS. H. LIMIT	2.000 Abs
NORMAL VOLUME	4 µL		
RERUN VOLUME	2 µL	REAGENT BLANK	Yes (0.000)
INCUBATION TIME	11.5 min	Cal. Low Limit	*
		FACTOR	#

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NAME	CÁLCIO CPC	PROZONE CHECK	No
ABBR. NAME	Ca	REF. MALE LOW	*
MODE	Endpoint	REF. MALE HIGH	*
WAVELENGTH	578 nm	REF. FEMALE LOW	*
UNITS	Mg/dL	REF. FEMALE HIGH	*
DECIMALS	2	REF. PED. LOW	*
LOW CONC.	0.00 mg/dL	REF. PED. HIGH	*
HIGH CONC.	20 mg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PANIC HIGH	*
REPEAT	1	CONTROL 1	*
NUMBER	1	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	CÁLCIO CPC		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	300 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	301 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	3 µL		
RERUN VOLUME	2 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	0µL	Cal. High Limit	*
RERUN VOLUME	0µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	4.5 min		
<b>MONO MODE</b>			
NAME	CÁLCIO CPC		
SAMPLE BLANK	No	LOW ABSORBANCE	0.000 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	300 µL	R. ABS. L LIMIT	1.100 Abs
RERUN VOLUME	301 µL	R. ABS. H. LIMIT	2.000 Abs
SAMPLE			
NORMAL VOLUME	3 µL	REAGENT BLANK	Yes (0.000)
RERUN VOLUME	2 µL	Cal. Low Limit	*
INCUBATION TIME	11.5 min	FACTOR	#

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NAME	CK MB	DELAY, MIN. TIME	285,293 sec
ABBR. NAME	CK MB		
MODE	Kinetic	PROZONE CHECK	No
WAVELENGTH	340 nm	REF. MALE LOW	*
UNITS	U/L	REF. MALE HIGH	*
DECIMALS	0	REF. FEMALE LOW	*
LOW CONC.	0 U/L	REF. FEMALE HIGH	*
HIGH CONC.	600 U/L	REF. PED. LOW	*
CALIBRATOR NAME	Utilizar Fator	REF. PED. HIGH	*
REPEAT	1	REF. PANIC LOW	*
NUMBER	1	REF. PANIC HIGH	*
CONCENTRATION	*	CONTROL 1	*
INTERVAL	0 days	CONTROL 2	*
CUT OFF	No	CONTROL 3	*
<b>DUAL MODE</b>		CORRELAT. FACTOR	1.000
NAME	CK MB	CORRELAT. OFFSET	0.00
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.600 Abs
NORMAL VOLUME	250 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	252 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	10 µL		
RERUN VOLUME	8 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	0 µL	Cal. High Limit	*
RERUN VOLUME	0 µL		
PREDILUTION	No	FACTOR	8254
INCUBATION TIME	289,230 sec		
<b>MONO MODE</b>		LOW ABSORBANCE	0.600 Abs
NAME	CK MB	HIGH ABSORBANCE	3.000 Abs
SAMPLE BLANK	No	R. ABS. L LIMIT	1.100 Abs
R1 BOTTLE	25 mL	R. ABS. H. LIMIT	2.000 Abs
NORMAL VOLUME	230 µL		
RERUN VOLUME	253 µL	REAGENT BLANK	Yes (0.000)
SAMPLE		Cal. Low Limit	*
NORMAL VOLUME	23 µL	FACTOR	8254
RERUN VOLUME	20 µL		

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NAME	CREATININA	INCUBATION TIME	32,188 sec
ABBR. NAME	CRE		
MODE	Twopoint	PROZONE CHECK	No
WAVELENGTH	505 nm	REF. MALE LOW	*
UNITS	mg/dL	REF. MALE HIGH	*
DECIMALS	1	REF. FEMALE LOW	*
LOW CONC.	0.0 mg/dL	REF. FEMALE HIGH	*
HIGH CONC.	20 mg/dL	REF. PED. LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PED. HIGH	*
REPEAT	1	REF. PANIC LOW	*
NUMBER	1	REF. PANIC HIGH	*
CONCENTRATION	*	CONTROL 1	*
INTERVAL	*	CONTROL 2	*
CUT OFF	No	CONTROL 3	*
<b>DUAL MODE</b>		CORRELAT. FACTOR	1.000
NAME	CREATININA	CORRELAT. OFFSET	0.00
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	250 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	270 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	25 µL		
RERUN VOLUME	5 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	0 µL	Cal. High Limit	*
RERUN VOLUME	0 µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	24,183 sec		
<b>MONO MODE</b>		LOW ABSORBANCE	0.000 Abs
NAME	CREATININA	HIGH ABSORBANCE	3.000 Abs
SAMPLE BLANK	No	R. ABS. L LIMIT	1.100 Abs
R1 BOTTLE	25 mL	R. ABS. H. LIMIT	2.000 Abs
NORMAL VOLUME	250 µL		
RERUN VOLUME	270 µL	REAGENT BLANK	Yes (0.000)
SAMPLE		Cal. Low Limit	*
NORMAL VOLUME	25 µL	FACTOR	#
RERUN VOLUME	5 µL		

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NAME	FERRO FERROZINE	PROZONE CHECK	No
ABBR. NAME	FER	REF. MALE LOW	*
MODE	End Point	REF. MALE HIGH	*
WAVELENGTH	546 nm	REF. FEMALE LOW	*
UNITS	µg/dL	REF. FEMALE HIGH	*
DECIMALS	0	REF. PED. LOW	*
LOW CONC.	0 µg/dL	REF. PED. HIGH	*
HIGH CONC.	500 µg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	Padrão FEZ	REF. PANIC HIGH	*
REPEAT	*	CONTROL 1	*
NUMBER	*	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	FERRO FERROZINE		
SAMPLE BLANK	Yes		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	200 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	220 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	40.0 µL		
RERUN VOLUME	20.0 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	50 µL	Cal. High Limit	*
RERUN VOLUME	30 µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	6.5		
<b>MONO MODE</b>			
NAME	FERRO FERROZINE		
SAMPLE BLANK	No	LOW ABSORBANCE	0.000 Abs
R1 BOTTLE	23 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	0 µL	R. ABS. L LIMIT	1.100 Abs
RERUN VOLUME	0 µL	R. ABS. H. LIMIT	2.000 Abs
SAMPLE			
NORMAL VOLUME	3 µL	REAGENT BLANK	Yes (0.000)
RERUN VOLUME	0 µL	Cal. Low Limit	*
INCUBATION TIME	0 min	FACTOR	#

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NAME	Ferro CRX	PROZONE CHECK	No
ABBR. NAME	Fe	REF. MALE LOW	*
MODE	Endpoint	REF. MALE HIGH	*
WAVELENGTH	620 nm	REF. FEMALE LOW	*
UNITS	µg/dL	REF. FEMALE HIGH	*
DECIMALS	0	REF. PED. LOW	*
LOW CONC.	0.00 µg/dL	REF. PED. HIGH	*
HIGH CONC.	500 µg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PANIC HIGH	*
REPEAT	1	CONTROL 1	*
NUMBER	1	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	Ferro CRX		
SAMPLE BLANK	No	LOW ABSORBANCE	0.000 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	240 µL	R. ABS. L LIMIT	- 0.100 Abs
RERUN VOLUME	245 µL	R. ABS. H. LIMIT	0.500 Abs
SAMPLE			
NORMAL VOLUME	10 µL	REAGENT BLANK	Yes (0.001)
RERUN VOLUME	5 µL	Cal. Low Limit	*
R2 BOTTLE	5 mL	Cal. High Limit	*
NORMAL VOLUME	0 µL		
RERUN VOLUME	0 µL	FACTOR	#
PREDILUTION	No		
INCUBATION TIME	4.5 min		
<b>MONO MODE</b>			
NAME	Ferro CRX		
SAMPLE BLANK	No	LOW ABSORBANCE	0.000 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	240 µL	R. ABS. L LIMIT	1.100 Abs
RERUN VOLUME	245 µL	R. ABS. H. LIMIT	2.000 Abs
SAMPLE			
NORMAL VOLUME	10 µL	REAGENT BLANK	Yes (0.000)
RERUN VOLUME	5 µL	Cal. Low Limit	*
INCUBATION TIME	11.5 min	FACTOR	#

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NAME	FÓSFORO UV	PROZONE CHECK	No
ABBR. NAME	FUV	REF. MALE LOW	*
MODE	Endpoint	REF. MALE HIGH	*
WAVELENGTH	340 nm	REF. FEMALE LOW	*
UNITS	mg/dL	REF. FEMALE HIGH	*
DECIMALS	2	REF. PED. LOW	*
LOW CONC.	0 mg/dL	REF. PED. HIGH	*
HIGH CONC.	15 mg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PANIC HIGH	*
REPEAT	1	CONTROL 1	*
NUMBER	1	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	FÓSFORO UV		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	260 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	262 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	4 µL		
RERUN VOLUME	2 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	0 µL	Cal. High Limit	*
RERUN VOLUME	0 µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	4.5 min		
<b>MONO MODE</b>			
NAME	FÓSFORO UV		
SAMPLE BLANK	No	LOW ABSORBANCE	0.000 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	260 µL	R. ABS. L LIMIT	1.100 Abs
RERUN VOLUME	262 µL	R. ABS. H. LIMIT	2.000 Abs
SAMPLE			
NORMAL VOLUME	4 µL	REAGENT BLANK	Yes (0.000)
RERUN VOLUME	2.0 µL	Cal. Low Limit	*
INCUBATION TIME	11.5 min	FACTOR	#

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NAME	FRUTOSAMINA	PROZONE CHECK	No
ABBR. NAME	FTS	REF. MALE LOW	*
MODE	Twopoint	REF. MALE HIGH	*
WAVELENGTH	546 nm	REF. FEMALE LOW	*
UNITS	mmol/L	REF. FEMALE HIGH	*
DECIMALS	2	REF. PED. LOW	*
LOW CONC.	0.0 mg/dL	REF. PED. HIGH	*
HIGH CONC.	7,0 mg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	Cal FTS	REF. PANIC HIGH	*
REPEAT NUMBER	1	CONTROL 1	*
CONCENTRATION	*	CONTROL 2	*
INTERVAL	*	CONTROL 3	*
CUT OFF	No	CORRELAT. FACTOR	1.000
<b>DUAL MODE</b>		CORRELAT. OFFSET	0.00
NAME	FRUTOSAMINA		
SAMPLE BLANK	No		
R1 BOTTLE	25 MI	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	250 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	270 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	1.500 Abs
NORMAL VOLUME	15 µL		
RERUN VOLUME	5 µL	REAGENT BLANK	Yes (0.000)
<b>MONO MODE</b>		Cal. Low Limit	*
NAME	FRUTOSAMINA	Cal. High Limit	*
SAMPLE BLANK	No		
R1 BOTTLE	25 ml	FACTOR	#
NORMAL VOLUME	225 µL		
RERUN VOLUME	225 µL		
SAMPLE			
NORMAL VOLUME	11 µL	LOW ABSORBANCE	0.000 Abs
RERUN VOLUME	25 µL	HIGH ABSORBANCE	3.000 Abs
PREDILUTION	No	R. ABS. H. LIMIT	2.000 Abs
SLOPE BLANK	No	REAGENT BLANK	Yes (0.000)
INCUBATION TIME	386, 519 sec	Cal. Low Limit	*
		FACTOR	#

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NAME	GAMA GT	PROZONE CHECK	No
ABBR. NAME	GGT	REF. MALE LOW	*
MODE	Kinetic	REF. MALE HIGH	*
WAVELENGTH	405 nm	REF. FEMALE LOW	*
UNITS	U/L	REF. FEMALE HIGH	*
DECIMALS	0	REF. PED. LOW	*
LOW CONC.	0 U/L	REF. PED. HIGH	*
HIGH CONC.	2500 U/L	REF. PANIC LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PANIC HIGH	*
REPEAT	1	CONTROL 1	*
NUMBER	1	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	GGT		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	240 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	250 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	20 µL		
RERUN VOLUME	10 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	0 µL	Cal. High Limit	*
RERUN VOLUME	0 µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	77,230 sec		
<b>MONO MODE</b>			
NAME	GGT		
SAMPLE BLANK	No	LOW ABSORBANCE	0.000 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	240 µL	R. ABS. L LIMIT	1.100 Abs
RERUN VOLUME	250 µL	R. ABS. H. LIMIT	2.000 Abs
SAMPLE			
NORMAL VOLUME	20 µL	REAGENT BLANK	Yes (0.000)
RERUN VOLUME	10 µL	Cal. Low Limit	*
DELAY, MIN. TIME	70,352 sec	FACTOR	#

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NAME	HDL COLESTEROL	PROZONE CHECK	No
ABBR. NAME	HDL	REF. MALE LOW	*
MODE	Endpoint	REF. MALE HIGH	*
WAVELENGTH	505 nm	REF. FEMALE LOW	*
UNITS	mg/dL	REF. FEMALE HIGH	*
DECIMALS	0	REF. PED. LOW	*
LOW CONC.	0 mg/dL	REF. PED. HIGH	*
HIGH CONC.	150 mg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	Padrão HDL	REF. PANIC HIGH	*
REPEAT	*	CONTROL 1	*
NUMBER	*	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	HDL COLESTEROL		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	250 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	301 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	25 µL		
RERUN VOLUME	2 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL		
PREDILUTION	No	Cal. Low Limit	*
INCUBATION TIME	4.5 min	Cal. High Limit	*
<b>MONO MODE</b>			
NAME	HDL COLESTEROL	FACTOR	#
SAMPLE BLANK	No	LOW ABSORBANCE	0.000 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	250 µL	R. ABS. L LIMIT	1.100 Abs
RERUN VOLUME	300 µL	R. ABS. H. LIMIT	2.000 Abs
SAMPLE			
NORMAL VOLUME	25 µL	REAGENT BLANK	Yes (0.000)
RERUN VOLUME	2 µL	Cal. Low Limit	*
INCUBATION TIME	11.5 min	FACTOR	#

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NAME	HDL DIRETO	PROZONE CHECK	No
ABBR. NAME	HDL DIRETO	REF. MALE LOW	*
MODE	Endpoint	REF. MALE HIGH	*
WAVELENGTH	620 nm	REF. FEMALE LOW	*
UNITS	mg/dL	REF. FEMALE HIGH	*
DECIMALS	0	REF. PED. LOW	*
LOW CONC.	0 mg/dL	REF. PED. HIGH	*
HIGH CONC.	180 mg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	Padrão HDL	REF. PANIC HIGH	*
REPEAT	1	CONTROL 1	*
NUMBER	1	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	HDL DIRETO		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	220 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	223 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	6 µL		
RERUN VOLUME	3 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	60 µL	Cal. High Limit	*
RERUN VOLUME	60 µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	6.5 min		

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NAME	LACTATO	PROZONE CHECK	No
ABBR. NAME	LAC	REF. MALE LOW	*
MODE	Endpoint	REF. MALE HIGH	*
WAVELENGTH	546 nm	REF. FEMALE LOW	*
UNITS	mg/dL	REF. FEMALE HIGH	*
DECIMALS	0	REF. PED. LOW	*
LOW CONC.	0 mg/dL	REF. PED. HIGH	*
HIGH CONC.	120 mg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PANIC HIGH	*
REPEAT	*	CONTROL 1	*
NUMBER	*	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	LACTATO		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	300 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	301 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	3 µL		
RERUN VOLUME	2 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	0 µL	Cal. High Limit	*
RERUN VOLUME	0 µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	11.5 min		
<b>MONO MODE</b>			
NAME	LACTATO		
SAMPLE BLANK	No	LOW ABSORBANCE	0.000 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	300 µL	R. ABS. L LIMIT	1.100 Abs
RERUN VOLUME	300 µL	R. ABS. H. LIMIT	2.000 Abs
SAMPLE			
NORMAL VOLUME	3 µL	REAGENT BLANK	Yes (0.000)
RERUN VOLUME	2 µL	Cal. Low Limit	*
INCUBATION TIME	11.5 min	FACTOR	#

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NAME	LDL DIRETO	PROZONE CHECK	No
ABBR. NAME	LDL DIRETO	REF. MALE LOW	*
MODE	Endpoint	REF. MALE HIGH	*
WAVELENGTH	578 nm	REF. FEMALE LOW	*
UNITS	mg/dL	REF. FEMALE HIGH	*
DECIMALS	0	REF. PED. LOW	*
LOW CONC.	0 mg/dL	REF. PED. HIGH	*
HIGH CONC.	400 mg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	Padrão LDL	REF. PANIC HIGH	*
REPEAT	1	CONTROL 1	*
NUMBER	1	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
DUAL MODE			
NAME	LDL DIRETO		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	220 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	223 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	6 µL		
RERUN VOLUME	3 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	60 µL	Cal. High Limit	*
RERUN VOLUME	60 µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	6.5 min		

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NAME	MAGNÉSIO	PROZONE CHECK	No
ABBR. NAME	MAG	REF. MALE LOW	*
MODE	Endpoint	REF. MALE HIGH	*
WAVELENGTH	505 nm	REF. FEMALE LOW	*
UNITS	mg/dL	REF. FEMALE HIGH	*
DECIMALS	0	REF. PED. LOW	*
LOW CONC.	0 mg/dL	REF. PED. HIGH	*
HIGH CONC.	4.5 mg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PANIC HIGH	*
REPEAT	*	CONTROL 1	*
NUMBER	*	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	MAGNÉSIO		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	300 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	301 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	4 µL		
RERUN VOLUME	2 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	0 µL	Cal. High Limit	*
RERUN VOLUME	0 µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	4.5 min		
<b>MONO MODE</b>			
NAME	MAGNÉSIO		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	300 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	300 µL	R. ABS. L LIMIT	1.100 Abs
SAMPLE		R. ABS. H. LIMIT	2.000 Abs
NORMAL VOLUME	4 µL		
RERUN VOLUME	2 µL	REAGENT BLANK	Yes (0.000)
INCUBATION TIME	11.5 min	Cal. Low Limit	*
		FACTOR	#

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NAME	PROTEÍNA URINÁRIA	PROZONE CHECK	No
ABBR. NAME	PRO-U	REF. MALE LOW	*
MODE	Endpoint	REF. MALE HIGH	*
WAVELENGTH	620 nm	REF. FEMALE LOW	*
UNITS	mg/L	REF. FEMALE HIGH	*
DECIMALS	0	REF. PED. LOW	*
LOW CONC.	0 mg/L	REF. PED. HIGH	*
HIGH CONC.	3000 mg/L	REF. PANIC LOW	*
CALIBRATOR NAME	Padrão Proteína Urinária	REF. PANIC HIGH	*
REPEAT	*	CONTROL 1	*
NUMBER	*	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	PROTEÍNA URINÁRIA		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	300 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	301 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	6 µL		
RERUN VOLUME	2 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	0 µL	Cal. High Limit	*
RERUN VOLUME	0 µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	11.5 min		
<b>MONO MODE</b>			
NAME	PROTEÍNA URINÁRIA		
SAMPLE BLANK	No	LOW ABSORBANCE	0.000 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	300 µL	R. ABS. L LIMIT	1.100 Abs
RERUN VOLUME	300 µL	R. ABS. H. LIMIT	2.000 Abs
SAMPLE			
NORMAL VOLUME	6 µL	REAGENT BLANK	Yes (0.000)
RERUN VOLUME	2 µL	Cal. Low Limit	*
INCUBATION TIME	11.5 min	FACTOR	#

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NAME	TRIGLICÉRIDES	PROZONE CHECK	No
ABBR. NAME	TGR	REF. MALE LOW	*
MODE	Endpoint	REF. MALE HIGH	*
WAVELENGTH	505 nm	REF. FEMALE LOW	*
UNITS	mg/dL	REF. FEMALE HIGH	*
DECIMALS	0	REF. PED. LOW	*
LOW CONC.	0 mg/dL	REF. PED. HIGH	*
HIGH CONC.	1000 mg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	AUTOCAL	REF. PANIC HIGH	*
REPEAT	*	CONTROL 1	*
NUMBER	*	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	TRIGLICÉRIDES		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	300 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	301 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	3 µL		
RERUN VOLUME	2 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	0 µL	Cal. High Limit	*
RERUN VOLUME	0 µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	11.5 min		
<b>MONO MODE</b>			
NAME	TRIGLICÉRIDES		
SAMPLE BLANK	No	LOW ABSORBANCE	0.000 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	300 µL	R. ABS. L LIMIT	1.100 Abs
RERUN VOLUME	300 µL	R. ABS. H. LIMIT	2.000 Abs
SAMPLE			
NORMAL VOLUME	3 µL	REAGENT BLANK	Yes (0.000)
RERUN VOLUME	2 µL	Cal. Low Limit	*
INCUBATION TIME	11.5 min	FACTOR	#

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NAME	URÉIA	PROZONE CHECK	No
ABBR. NAME	UUU	REF. MALE LOW	*
MODE	Twopoint	REF. MALE HIGH	*
WAVELENGTH	340 nm	REF. FEMALE LOW	*
UNITS	mg/dL	REF. FEMALE HIGH	*
DECIMALS	1	REF. PED. LOW	*
LOW CONC.	0.0 mg/dL	REF. PED. HIGH	*
HIGH CONC.	250 mg/dL	REF. PANIC LOW	*
CALIBRATOR NAME	AUTOCAL H	REF. PANIC HIGH	*
REPEAT	1	CONTROL 1	*
NUMBER	1	CONTROL 2	*
CONCENTRATION	*	CONTROL 3	*
INTERVAL	*	CORRELAT. FACTOR	1.000
CUT OFF	No	CORRELAT. OFFSET	0.00
<b>DUAL MODE</b>			
NAME	URÉIA		
SAMPLE BLANK	No		
R1 BOTTLE	25 mL	LOW ABSORBANCE	0.000 Abs
NORMAL VOLUME	300 µL	HIGH ABSORBANCE	3.000 Abs
RERUN VOLUME	255 µL	R. ABS. L LIMIT	- 0.100 Abs
SAMPLE		R. ABS. H. LIMIT	0.500 Abs
NORMAL VOLUME	4 µL		
RERUN VOLUME	2 µL	REAGENT BLANK	Yes (0.001)
R2 BOTTLE	5 mL	Cal. Low Limit	*
NORMAL VOLUME	0 µL	Cal. High Limit	*
RERUN VOLUME	0 µL		
PREDILUTION	No	FACTOR	#
INCUBATION TIME	24,120 sec		
<b>MONO MODE</b>			
NAME	URÉIA		
SAMPLE BLANK	No	LOW ABSORBANCE	0.000 Abs
R1 BOTTLE	25 mL	HIGH ABSORBANCE	3.000 Abs
NORMAL VOLUME	300 µL	R. ABS. L LIMIT	1.100 Abs
RERUN VOLUME	300 µL	R. ABS. H. LIMIT	2.000 Abs
SAMPLE			
NORMAL VOLUME	4 µL	REAGENT BLANK	Yes (0.000)
RERUN VOLUME	2 µL	Cal. Low Limit	*
INCUBATION TIME	32,120 sec	FACTOR	#

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